

Development and Validation of an HPLC Method for Simultaneous Determination of Nortriptyline and Pregabalin in Pharmaceutical Dosage Form

Karuru Anil^{*1}, R. Kiran Jyothi², M. Mahesh³

^{*1}Student, Department of Pharmaceutical analysis, Oil Technological & Pharmaceutical Research Institute, JNTU, Anantapur, Andhra Pradesh, India.

^{2,3} Department of Pharmaceutical analysis, Oil Technological & Pharmaceutical Research Institute, JNTU, Anantapur, Andhra Pradesh, India.

Submitted: 15-08-2023

Accepted: 25-08-2023

ABSTRACT: The assay of Pregabalin and Nortriptyline was performed with tablets and the % assay was found to be 99.58 and 99.61 which shows that the method is useful for routine analysis.

The linearity of Pregabalin and Nortriptyline was found to be linear with a correlation coefficient of 0.999 and 0.999, which shows that the method is capable of producing good sensitivity.

The acceptance criteria of precision is RSD should be not more than 2.0% and the method show precision 0.33 and 0.27 for Pregabalin and Nortriptyline which shows that the method is precise.

The acceptance criteria of intermediate precision is RSD should be not more than 2.0% and the method show precision 0.47 and 0.30 for Pregabalin and Nortriptyline which shows that the method is repeatable when performed in different days also.

The accuracy limit is the percentage recovery should be in the range of 97.0% - 103.0%. The total recovery was found to be 99.09% and 100.73% for Pregabalin and Nortriptyline. The validation of developed method shows that the accuracy is well within the limit, which shows that the method is capable of showing good accuracy and reproducibility.


KEYWORDS: Pregabalin, Nortriptyline, RSD, accuracy

I. INTRODUCTION

PREGABALIN:

Pregabalin is an anticonvulsant drug used to treat neuropathic pain conditions and fibromyalgia, and for the treatment of partial onset seizures in combination with other anticonvulsants.

Chemical structure:



IUPACNAME: (3S)-3-(aminomethyl)-5-methylhexanoic acid

Chemical Formula: C₈H₁₇NO₂

Molecular weight: 159.22


Melting point: 176 - 178°C

Pka: 4.8

2. NORTRIPTYLINE:

Description: Nortriptyline is a tricyclic antidepressant used in the treatment of depression

Chemical Structure:



Formula : C₁₉H₂₁N

Molecular weight : 263.37

Freely soluble : Methanol

PHARMACOLOGY:

Pharmacodynamics:

Nortriptyline exerts antidepressant effects likely by inhibiting the reuptake of serotonin and norepinephrine at neuronal cell membranes. It also exerts antimuscarinic effects through its actions on the acetylcholine receptor.

Mechanism of action:

Though prescribing information does not identify a specific mechanism of action for nortriptyline¹⁸, is believed that nortriptyline either inhibits the reuptake of the neurotransmitter serotonin at the neuronal membrane or acts at the level of the beta-adrenergic receptors. It displays a more selective reuptake inhibition for noradrenaline, which may explain increased symptom improvement after nortriptyline therapy.¹⁶ Tricyclic antidepressants do not inhibit monoamine oxidase nor do they affect dopamine reuptake.¹⁸ As with other tricyclics, nortriptyline displays affinity for other receptors including mACh receptors, histamine receptors, 5-HT receptors, in addition to other receptors

Brand Names:



II. METHODOLOGY

HPLC METHOD DEVELOPMENT:

Wave length selection:

UV spectrum of 10 µg / ml Pregabalin and Nortriptyline in diluents (mobile phase composition) was recorded by scanning in the range of 200nm to 400nm. From the UV spectrum wavelength selected as 227. At this wavelength both the drugs show good absorbance.

Optimization of Column: Waters Symmetry C₁₈ (4.6 x 250mm, 5.0µm) was found to be ideal as it gave good peak shape and resolution at 1.2 ml/min flow.

OPTIMIZED CHROMATOGRAPHIC CONDITIONS:

Instrument used : High performance liquid chromatography equipped with Auto Sampler and PDA detector

Temperature: Ambient

Column : Waters Symmetry C₁₈ (4.6 x 250mm, 5.0µm)

Buffer: 0.1% ortho phosphoric acid buffer

Mobile phase: 25% buffer: 75% Methanol

Flow rate: 1.0 ml per min

Wavelength: 227 nm

Injection volume : 20 µl

Run time: 10min.

Optimized chromatogram is shown in the **fig 2**. System suitability parameters are shown in **fig 2** and the results are shown in **Table 6**.

PREPARATION OF BUFFER AND MOBILE PHASE:

Preparation of 0.1% Ortho phosphoric acid buffer:

Pipetted 1 ml of ortho phosphoric acid in 100 ml HPLC water.

Preparation of mobile phase:

Mix a mixture of above buffer 250 ml (25%) and 750 ml Methanol HPLC (75%) and degas in ultrasonic water bath for 5 minutes. Filter through 045 µ filter under vacuum filtration.

Diluent Preparation:

Use the Mobile phase as Diluents.

VALIDATION PARAMETERS:

1. ASSAY:

Standard Solution Preparation:

Accurately weigh and transfer 75mg of Pregabalin & 10mg of Nortriptyline working standard into a 100ml clean dry volumetric flask add Diluents and sonicate to dissolve it completely and make volume up to the mark with the same solvent. (Stock solution)

Further pipette 1.5 ml of Pregabalin & Nortriptyline of the above stock solution into a 10ml volumetric flask and dilute up to the mark with Diluents.

Sample Solution Preparation:

Accurately weigh and transfer equivalent to 75mg of Pregabalin & 10mg Nortriptyline equivalent weight of the sample into a 100ml clean dry volumetric flask add about 70ml of Diluents and sonicate to dissolve it completely and make volume up to the mark with the same solvent. (Stock solution).

Further pipette 1.5 ml of Pregabalin & Nortriptyline of the above stock solution into a 10ml volumetric flask and dilute up to the mark with Diluents.

Procedure: Inject 20 μ L of the standard, sample into the chromatographic system and measure the areas for the Pregabalin & Nortriptyline peaks and calculate the %Assay by using the formulae. The chromatograms were recorded as show in **Fig 3, 4** and results are shown in **Table- 7**.

2. LINEARITY:

Preparation of stock solution:

Accurately weigh and transfer 75mg of Pregabalin & 10mg of Nortriptyline working standard into a 100ml clean dry volumetric flask add about 70mL of Diluents and sonicate to dissolve it completely and make volume up to the mark with the same solvent. (Stock solution).

Preparation of Level – I (37.5 ppm & 5 ppm of Pregabalin & Nortriptyline):

0.5 ml of stock solution has taken in 10ml of volumetric flask dilute up to the mark with Diluents.

Preparation of Level – II (75 ppm & 10 ppm of Pregabalin & Nortriptyline):

1.0 ml of stock solution has taken in 10ml of volumetric flask dilute up to the mark with Diluents.

Preparation of Level – III (112.5 ppm & 15 ppm of Pregabalin & Nortriptyline):

1.5 ml of stock solution has taken in 10ml of volumetric flask dilute up to the mark with Diluents.

Preparation of Level – IV (150 ppm & 20ppm of Pregabalin & Nortriptyline):

2.0 ml of stock solution has taken in 10ml of volumetric flask dilute up to the mark with Diluents.

Preparation of Level – V (187.5 ppm & 25 ppm of Pregabalin & Nortriptyline):

5ml of stock solution has taken in 10ml of volumetric flask dilute up to the mark with Diluents.

Procedure: Inject each level into the chromatographic system and measure the peak area. Plot a graph of peak area versus concentration (on X-axis concentration and on Y-axis Peak area) and calculate the correlation coefficient.

The chromatograms were recorded as show in **Fig 5 – 9**, calibrationgraphs were shown in **Fig 10,11** and results are shown in **Table-8, 9**.

3. PRECISION:

Preparation of stock Solution:

Accurately weigh and transfer 75mg of Pregabalin & 10mg of Nortriptyline working standard into a 100ml clean dry volumetric flask add

about 70ml of Diluents and sonicate to dissolve it completely and make volume up to the mark with the same solvent. (Stock solution)

Further pipette 1.5 ml of Pregabalin & Nortriptyline the above stock solution into a 10ml volumetric flask and dilute up to the mark with Diluents.

Procedure: The standard solution was injected for six times and measured the area for all five injections in HPLC. The %RSD for the area of five replicate injections was found to be within the specified limits.

The sample solutions of Pregabalin and Nortriptyline chromatograms were recorded as shown in **Fig 12- 16**. The mean and percentage relative standard deviation were calculated from the peak areas and shown in the **Table 10 and 11**.

4. Intermediate Precision/Ruggedness: To evaluate the intermediate precision (also known as Ruggedness) of the method, Precision was performed on different day within the laboratory.

Preparation of stock solution:

Accurately weigh and transfer 75mg of Pregabalin & 10mg of Nortriptyline working standard into a 100ml clean dry volumetric flask add about 70ml of Diluents and sonicate to dissolve it completely and make volume up to the mark with the same solvent. (Stock solution)

Further pipette 1.5ml of Pregabalin & Nortriptyline the above stock solution into a 10ml volumetric flask and dilute up to the mark with Diluents.

Procedure: The standard solution was injected for Six times and measured the area for all Five injections in HPLC. The %RSD for the area of five replicate injections was found to be within the specified limits.

Chromatograms were recorded as shown in **Fig 17- 22** and results are shown in **Table 12 and 13**.

5. ACCURACY:

For accuracy determination, three different concentrations were prepared separately i.e. 50%, 100% and 150% for the analyte and chromatograms are recorded for the same.

Preparation of Standard stock solution:

Accurately weigh and transfer 75mg of Pregabalin & 10mg of Nortriptyline working standard into a 100ml clean dry volumetric flask add about 70ml of Diluents and sonicate to dissolve it completely and make volume up to the mark with the same solvent. (Stock solution)

Further pipette 1.5 ml of Pregabalin & Nortriptyline the above stock solution into a 10ml volumetric flask and dilute up to the mark with Diluents.

Preparation Sample solutions:

For preparation of 50% solution (With respect to target Assay concentration):

Accurately weigh and transfer 37.5mg of Pregabalin and 5mg of Nortriptyline working standard into a 100ml clean dry volumetric flask add about 70 ml of Diluents and sonicate to dissolve it completely and make volume up to the mark with the same solvent. (Stock solution). Further pipette 1.5 ml of Pregabalin & Nortriptyline the above stock solution into a 10ml volumetric flask and dilute up to the mark with Diluents

For preparation of 100% solution (With respect to target Assay concentration):

Accurately weigh and transfer 75mg of Pregabalin and 10mg of Nortriptyline working standard into a 100ml clean dry volumetric flask add about 70ml of Diluents and sonicate to

dissolve it completely and make volume up to the mark with the same solvent. (Stock solution).

Further pipette 1.5 ml of Pregabalin & Nortriptyline the above stock solution into a 10ml volumetric flask and dilute up to the mark with Diluents.

For preparation of 150% solution (With respect to target Assay concentration):

Accurately weigh and transfer 112.5mg of Pregabalin and 15mg of Nortriptyline equivalent weight of tablet powder into a 100ml clean dry volumetric flask add about 7ml of Diluents and sonicate to dissolve it completely and make volume up to the mark with the same solvent. (Stock solution). Further pipette 1.5 ml of Pregabalin & Nortriptyline the above stock solution into a 10ml volumetric flask and dilute up to the mark with Diluents.

TABLE 4: INSTRUMENTS USED

SL. No	Instrument	Model
1	HPLC	WATERS, software: Empower, 2695 separation module.2487 UV detector.
2	UV/VIS spectrophotometer	LABINDIA UV 3000 ⁺
3	pH meter	Adwa – AD 1020
4	Weighing machine	Afcoset ER-200A
5	Pipettes and Burettes	Borosil
6	Beakers	Borosil

TABLE 5: CHEMICALS USED

SL.No	Chemical	Brand
1	Pregabalin	Supplied by Kenna Helath care
2	Nortriptyline	Supplied by Kenna Helath care
3	Ortho phosphoric acid	FINAR chemical LTD
4	Water and Methanol for HPLC	Standard solutions Ltd
5	Acetonitrile for HPLC	Standard solutions Ltd
6	HCl, H ₂ O ₂ , NaOH	MERCK

III. RESULTS AND DISCUSSION




FIGURE 2: CHROMATOGRAM FOR SYSTEM SUITABILITY

TABLE 6: RESULTS OF SYSTEM SUITABILITY PARAMETERS

S.No	Name	RT(min)	Area (μV sec)	Height (μV)	USP resolution	USP tailing	USP count plate
1	Pregabalin	2.401	2094603	196622		1.71	2947.68
2	Nortriptyline	3.374	3694090	286174	4.38	1.61	3826.77

Acceptance criteria:

- Resolution between two drugs must be not less than 2.
- Theoretical plates must be not less than 2000.
- Tailing factor must be not more than 2.
- It was found from above data that all the system suitability parameters for developed method were within the limit.

VALIDATION PARAMETERS:

1. ASSAY:

Standard and sample solution injected as described under experimental work. The corresponding chromatograms and results are shown below.




FIGURE 3: CHROMATOGRAM FOR STANDARD




FIGURE 4: CHROMATOGRAM FOR SAMPLE

TABLE 7: RESULTS OF ASSAY FOR PREGABALIN AND NORTRIPTYLINE

	LABEL CLAIM (MG)	% ASSAY
PREGABALIN	75	99.58
NORTRIPTYLINE	10	99.61

2. LINEARITY:

THE LINEARITY RANGE WAS FOUND TO LIE FROM 1µG/ML TO 5µG/ML OF

PREGABALIN, 7.5µG/ML TO 37.5µG/ML OF NORTRIPTYLINE AND CHROMATOGRAMS ARE SHOWN BELOW.




FIGURE 5: CHROMATOGRAM FOR LINEARITY-1




FIGURE 6: CHROMATOGRAM FOR LINEARITY-2




FIGURE 7: CHROMATOGRAM FOR LINEARITY-3




FIGURE 8: CHROMATOGRAM FOR LINEARITY-4




FIGURE 10: CALIBRATION GRAPH FOR PREGABALIN




FIGURE 11: CALIBRATION GRAPH FOR NORTRIPTYLINE




FIGURE 32: CHROMATOGRAM OF PREGABALIN, NORTRIPTYLINE SHOWING LOD

TABLE 16: RESULTS OF LOD

Drug name	Baseline noise(μV)	Signal obtained (μV)	S/N ratio
Pregabalin	55	164	2.98
Nortriptyline	55	162	2.95




FIGURE 33: CHROMATOGRAM OF PREGABALIN, NORTRIPTYLINE SHOWING LOQ

TABLE 17: RESULTS OF LOQ

Drug name	Baseline noise(μV)	Signal obtained (μV)	S/N ratio
Pregabalin	55	548	9.96
Nortriptyline	55	549	9.98

IV. CONCLUSION

The accuracy limit is the percentage recovery should be in the range of 97.0% - 103.0%. The total recovery was found to be 99.09% and 100.73% for Pregabalin and Nortriptyline. The validation of developed method shows that the accuracy is well within the limit, which shows that the method is capable of showing good accuracy and reproducibility.

The acceptance criteria for LOD and LOQ is 3 and 10. The LOD and LOQ for Pregabalin was found to be 2.98 and 9.96 and LOD and LOQ for Nortriptyline was found to be 2.95 and 9.98.

The robustness limit for mobile phase variation and flow rate variation are well within the limit, which shows that the method is having good system suitability and precision under given set of conditions. The developed method can be used in the routine drug analysis of pharmaceutical formulations.

REFERENCES

- [1]. Gurudeep Chatwal and Sham Anand, Instrumental methods of chemical analysis. Himalaya publishers, 7th edition, 1992, Pg. no 2.624-2.639.
- [2]. Skoog et al., Principles of Instrumental Analysis. Barkhanath publishers, 8th edition, Pg. no. 973-995.
- [3]. Hobart.H.Willard et al., Instrumental methods of analysis, CBS Publications and Distributors, New Delhi, 1st edition ,1986, Pg. no 529-563.
- [4]. Sethi P.D., Quantitative analysis of Drugs & Pharmaceuticals. CBS publishers and distributors, New Delhi, 3rd edition, 2001 Pg. no 1-120.
- [5]. D.Helmeste et al., J.Chromatogr. Published in 1997, Pg. no 195-201.
- [6]. Interna ICH of technical requirements for the registration of pharmaceuticals for human use, validation of analytical parameters; methodology adopted in 1996, Geneva.
- [7]. ICH Guidelines Q2B, Validation of Analytical Procedure: Definitions, published in March 1996, Geneva, Switzerland.
- [8]. Singh GN, Gupta RP. Stability of Pharmaceuticals. Dept. of Pharmaceutics, Institute of Technology, B.H.U., Varanasi U.P. Eastern Pharmacist Aug1987; 85-89.
- [9]. Teresa I, Lucas, Bishara RH, Robert HS. A Stability Program for the Distribution of Drug Products. Pharmaceutical Technology Resource Guide 2004; 86-9.
- [10]. Gennaro AR. Remington: The Science and Practice of Pharmacy. Philadelphia, 12th Edn 2000; 986-90.
- [11]. Singh S, Bakshi M. Development of validated stability-indicating assay methods critical review. J Pharm Biomed Anal 2002; 28: 1011-40.
- [12]. ICH Harmonised Tripartite Guideline Stability Testing of New Drug Substances and Products Q1A (R2) 2003; 2: 1-24.
- [13]. Haritha Potluri 1,* , Sreenivasa Rao Battula 2, Sunandamma Yeturu Validated Stability Indicating Rp-Hplc Method For Simultaneous Determination Of Nortriptyline And Pregabalin In Bulk And Combined Dosage Formulations, J. Chil. Chem. Soc. Vol.62 No.2 Concepción Jun. 2017.
- [14]. Patel, Rajesh & Dholakiya, Sandip & Golwala, Dharmesh & Gohil, Naitik & Vaidya, Santosh. (2020). A Novel Spectrophotometric and RP-HPLC methods for Determination of Nortriptyline hydrochloride and Pregabalin in Tablets. 9. 1-9.
- [15]. Mevada, N.A., Patel, B.R., Patel, J.G., Vegad, K.L., Patel, V.S. (2017). Stability Indicating RP-HPLC Method Development and Validation for Simultaneous Estimation of Pregabalin and Nortriptyline in Tablet. 6(2), 1-7.
- [16]. Gadarajan TVS, Arunkumar N, Aiyalu R, Ponnilaravasan I, Tamilselvi N. Development and Validation of RP-HPLC Method for the Simultaneous Estimation of Gabapentin and Nortriptyline Hydrochloride in Bulk and Tablet Dosage Form. J Young Pharm. 2022;14(2):187-91.
- [17]. J.Anil Mohan, B.Rajkumar, T.Bhavya, A.Ashok Kumar, Rp-Hplc Method Development And Validation For The Simultaneous Quantitative Estimation Of Pregabalin, Mecobalamin And Alpha Lipoic Acid In Capsules, Int J Pharm Pharm Sci, Vol 6, Issue 1, 270-277.